EN	Trima Accel® Disposable Tubing Sets Instructions for Use Rx Only
BG	Набори тръби Trima Accel® за еднократна употреба Инструкции за употреба
CS	Jednorázové sety Trima Accel® Návod k použití
DA	Trima Accel® engangsslangesæt Brugervejledning
DE	Trima Accel®-Einwegschlauchsets Gebrauchsanweisung
EL	Αναλώσιμα σετ σωλήνων Trima Accel® Οδηγίες χρήσης
ES	Equipos de líneas desechables de Trima Accel® Instrucciones de uso
FI	Kertakäyttöiset Trima Accel® -letkusarjat Käyttöohjeet
FR	Kits de collecte à usage unique Trima Accel® Mode d'emploi
HR	Set cjevčica Trima Accel® za jednokratnu upotrebu Upute za upotrebu
HU	Trima Accel® egyszer használatos csőszerszélékek Használati utasítás
IT	Circuiti monouso Trima Accel® Istruzioni per l'uso
LT	„Trima Accel®“ vienkartiniai vamzdelių rinkiniai Naudojimo instrukcijos
NL	Trima Accel® disposable slangensets Gebruiksaanwijzing
NO	Trima Accel® engangsslangesett Indikasjoner for bruk
PL	Jednorazowe zestawy drenów Trima Accel® Instrukcja obsługi
PT	Conjuntos Descartáveis Trima Accel® Instruções de Uso
RO	Seturile de tubulatură de unică folosință Trima Accel® Instrucțiuni de utilizare
RU	Одноразовые наборы магистралей Trima Accel® Инструкция по эксплуатации
SV	Trima Accel® engångsslangset Bruksanvisning
TR	Trima Accel Tek Kullanımlık Setler Kullanım Yönergeleri
UK	Одноразові комплекти магістралей Trima Accel® Інструкції із застосування

82310

EN

Trima Accel® LRS® Platelet + Sampler + Auto PAS, MultiPlasma Set

BG

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Trima Accel® disposable tubing sets provide the interface between the donor and the Trima Accel Automated Blood Collection System. Depending on the tubing set used, platelets, plasma and/or red blood cells may be collected.

- Adequate studies have not been performed to evaluate the quality of gamma irradiation or freezing of ACD-A/AS-3 red blood cell (RBC) products collected with gravity drain leukoreduction process (TLR filter) on the Trima Accel system.
- Studies have not been performed to support gamma irradiation or freezing of ACD-A/AS-3 RBC products collected with an integrated in-line RBC leukoreduction filter(s) (Auto RBC filter) on the Trima Accel system.

The Trima Accel tubing sets are intended to route and collect blood components from the access site through the Trima Accel Automated Blood Collection System.

Please refer to the Trima Accel Automated Blood Collection system operator's manual for a complete listing of contraindications, warnings, cautions, adverse effects, operating instructions, clinical benefits, and performance characteristics for the Trima Accel system.

Contraindications

There are no known contraindications for use of the Trima Accel system, except those associated with the infusion of solutions as required by the procedure and those associated with all types of automated blood component separators.

Adverse Events

Any adverse event or serious incident that occurs in relation to this device should be reported to the local competent authority and to Terumo BCT. Customers in the U.S. should contact Terumo BCT Customer Support. Customers outside the U.S. should contact their local Terumo BCT representative.

Warnings

1. If the integrity of the tubing set is compromised for any reason, the tubing set is no longer functionally closed, and the product should not be stored beyond 24 hours.
2. Check that all pump headers have been properly loaded before beginning a procedure. Visually inspect each pump and carefully run your finger around the edge of each pump.
3. Terumo BCT, Inc. products bearing the Hazardous Substances symbol contain the following substance defined as carcinogenic, mutagenic or reprotoxic (CMR), category 1A or 1B, and/or a substance having endocrine-disrupting properties in a concentration above 0.1% weight by weight:
 - Bis(2-ethylhexyl) phthalate (DEHP), which is reprotoxic 1B (or CMR 1B) and has endocrine-disrupting properties; CAS No.: 117-81-7; EC No.: 204-211-0
 - Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

For more information, visit the European Chemicals Agency website at <https://echa.europa.eu/home>.

Whole blood donors are not exposed to DEHP. The potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk for potential harmful effects of exposure to DEHP. However, regulatory bodies have noted

that the benefit of doing a needed procedure may justify the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.

The needle contains cobalt in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

4. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

5. The blood and fluid pathways of the tubing set are sterilized with ethylene oxide and are nonpyrogenic. The packaging cover of the tubing set serves to protect the set and keep contents in place during shipping. Under any of the following conditions, DO NOT use the set:
- The end caps are not in place
 - The needle cap plug is not in place
 - Severe kinks in the tubing are apparent
 - The tubing set is incorrectly assembled
 - The tubing set is damaged
 - The clamps are closed on the pre-attached filter tubing assembly

Cautions

1. Each operator should be thoroughly familiar with the Trima Accel system operator's manual. All procedures should be performed by qualified medical personnel under the supervision of a physician.
2. Use aseptic technique throughout all procedures to ensure donor safety and product quality.

3. When loading the channel and the LRS chamber into the centrifuge filler, be careful not to stretch the tubing, as leaks could develop.
4. Be careful to press on tubing only with your fingers when loading the channel. To avoid accidental puncture, never use any sharp object.
5. When Trima Accel Auto RBC Sets are used, RBC products collected from certain donors may have the potential for ineffective filtration and leukoreduction.
6. When using MultiPlasma tubing sets, ensure that at least two plasma line clamps are open for procedures with plasma volumes greater than 600 mL.

Procedures for Use

Refer to the Trima Accel Automated Blood Collection System operator's manual for detailed information regarding the following procedures:

- Setting up the disposable tubing set
- Performing a collection procedure
- Removing products and the disposable tubing set

Supplemental Instructions for Connecting Anticoagulant (AC)

When you are prompted to connect AC, connect the luer on the AC line (orange) to the AC bag. Break the frangible connector on the AC bag to allow the AC to flow. Squeeze the drip chamber in order to fill it approximately halfway. Follow the system prompts to load the AC line into the AC sensor, and continue according to the Trima Accel system operator's manual.

Note: When you break frangible connectors, bend them in both directions to ensure that you break them completely. Failure to do so may result in restricted flow.

Clinical Benefits

Terumo BCT does not make any claims regarding benefits for the patients treated with blood products produced using the Trima Accel system. The Trima Accel system is an automated blood cell separator intended for use in the collection of blood components for transfusion into patients. Blood and blood component transfusion are widely used in locations throughout the world to provide lifesaving care to patients. Trima Accel facilities are making the benefits of blood donation available to patients who may need them for blood transfusion procedures. Collecting blood components via apheresis offers additional benefits such as targeted collections of specific components while returning other whole blood components back to the donor.

Target Donor Population

The target donor population is based on the type of apheresis procedure performed. Donor populations include individuals who meet local requirements for donation of the applicable blood component or components.

Storage Conditions

0 °C to 30 °C (32 °F to 86 °F)

Specifications

Tubing Set Components

Characteristic		Performance	Conditions
Primary Tubing Material		PVC plasticized with DEHP	
Access Needle Size	Gauge	17 (for 82310, 82410, 82420, 82460, 82470, 82616, and 82656)	
		18 (for 82520)	
	Cannula Length	24.3 mm (0.96 in)	

Tubing Specifications (non-sterile, nominal* dimensions)

Item	Wall Thickness	Outer Diameter	Inner Diameter
Blue Sleeve Tubing on Sample Bag	0.94 mm (0.037 in)	4.75 mm (0.187 in)	2.87 mm (0.113 in)
Tubing to Plasma Bag(s)	0.56 mm (0.022 in)	4.06 mm (0.160 in)	2.94 mm (0.116 in)
Tubing to Platelet Bags (82410, 82420, 82460, 82470, 82616, and 82656 only)	0.56 mm (0.022 in)	4.06 mm (0.160 in)	2.94 mm (0.116 in)
Clear Tubing to RBC Bags (82410, 82420, and 82460 only)	0.56 mm (0.022 in)	4.06 mm (0.160 in)	2.94 mm (0.116 in)
Numbered Tubing to RBC Bags (82410, 82420, and 82460 only)	0.61 mm (0.024 in)	4.19 mm (0.165 in)	3.00 mm (0.118 in)
Auto RBC Filter Kit Tubing and Auto RBC Additive Solution Line: Tubing from Frangible (Breakable) Connector to Cassette (82470, 82520, 82616, and 82656 only)	0.56 mm (0.022 in)	4.06 mm (0.160 in)	2.94 mm (0.116 in)
Numbered Tubing to Auto RBC Bags (82470, 82520, 82616, and 82656 only)	0.61 mm (0.024 in)	4.19 mm (0.165 in)	3.00 mm (0.118 in)
Access Needle Tubing	1.00 mm (0.040 in)	5.00 mm (0.197 in)	3.00 mm (0.118 in)
Blood Component Sampling Set: Tubing from Y Connector to Frangible (Breakable) Connector (82460 and 82616 only)	0.66 mm (0.026 in)	4.19 mm (0.165 in)	2.87 mm (0.113 in)
Auto PAS Additive Solution Line: Tubing from Sterile Barrier Filter to Check Valve (82310, 82410, 82420, 82460, and 82470 only)	0.94 mm (0.037 in)	4.75 mm (0.187 in)	2.87 mm (0.113 in)

***Note:** The tubing dimensions listed on this document are nominal values and are intended for use in determining compatibility with various types of laboratory equipment. The nominal values are specified target dimensions; however, because they are based on non-sterile tubing measurements, and due to variations in the manufacturing process, the actual dimensions may be slightly different.

Return of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department.

IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.

Please contact your local representative for information regarding returned goods and product complaints.

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©2023 Terumo BCT, Inc.
Terumo BCT, Inc. (“Terumo Blood and Cell Technologies”, “Terumo BCT”)
Terumo Russia LLC